

This assay sheet is also available on :

[www.medonic.se](http://www.medonic.se)[www.exigo-vet.com](http://www.exigo-vet.com)

LOT: 11701-52  
EXP: 2017-06-26  
Open Vial Stability: 16 days

<b>Medonic CA530, CA620</b>		<b>Mean</b>	<b>± Range</b>	<b>Expected Range</b>	
<b>WBC</b>	10 <sup>9</sup> /L	8.4	0.6	7.8	9.0
<b>Hct</b>	%	31.8	2.5	29.3	34.3
<b>MCV</b>	fL	80.0	5.0	75.0	85.0
<b>RBC</b>	10 <sup>12</sup> /L	3.98	0.18	3.80	4.16
<b>Hgb</b>	g/dL	11.7	0.4	11.3	12.1
<b>MCH</b>	pg	29.4	1.5	27.9	30.9
<b>MCHC</b>	g/dL	36.8	2.2	34.6	39.0
<b>PLT</b>	10 <sup>9</sup> /L	238	30	208	268

<b>EXIGO</b>		<b>Mean</b>	<b>± Range</b>	<b>Expected Range</b>	
<b>WBC</b>	10 <sup>9</sup> /L	8.0	0.6	7.4	8.6
<b>Hct</b>	%	32.5	2.5	30.0	35.0
<b>MCV</b>	fL	82.6	5.0	77.6	87.6
<b>RBC</b>	10 <sup>12</sup> /L	3.93	0.18	3.75	4.11
<b>Hgb</b>	g/dL	11.7	0.4	11.3	12.1
<b>MCH</b>	pg	29.8	1.5	28.3	31.3
<b>MCHC</b>	g/dL	36.0	2.2	33.8	38.2
<b>PLT</b>	10 <sup>9</sup> /L	211	30	181	241

Manufactured for:  
Boule Medical AB  
Domnarvsgatan 4  
SE-163 53 Spånga, Sweden



1504056

**For EXIGO/EXIGO EOS with barcode reader**

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"

2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo  
Lot: 11701-52  
Prod.: Vet-Con Norm  
ID: 1170152+  
Exp.: 26-Jun-2017

