

LOT: 11710-52  
EXP: 2018-03-26  
Open Vial Stability: 16 days

Medonic CA530, CA620		Mean	± Range	Expected Range	
WBC	10 <sup>9</sup> /L	8,3	0,6	7,7	8,9
Hct	%	28,9	2,5	26,4	31,4
MCV	fL	74,1	5,0	69,1	79,1
RBC	10 <sup>12</sup> /L	3,90	0,18	3,72	4,08
Hgb	g/dL	10,8	0,4	10,4	11,2
MCH	pg	27,7	1,5	26,2	29,2
MCHC	g/dL	37,4	2,2	35,2	39,6
PLT	10 <sup>9</sup> /L	247	30	217	277

EXIGO		Mean	± Range	Expected Range	
WBC	10 <sup>9</sup> /L	7,8	0,6	7,2	8,4
Hct	%	29,4	2,5	26,9	31,9
MCV	fL	76,5	5,0	71,5	81,5
RBC	10 <sup>12</sup> /L	3,84	0,18	3,66	4,02
Hgb	g/dL	10,9	0,4	10,5	11,3
MCH	pg	28,4	1,5	26,9	29,9
MCHC	g/dL	37,1	2,2	34,9	39,3
PLT	10 <sup>9</sup> /L	193	30	163	223

This assay sheet is also available on :

[www.medonic.se](http://www.medonic.se)  
[www.exigo-vet.com](http://www.exigo-vet.com)

## For EXIGO/EXIGO EOS with barcode reader

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Instr.: Exigo  
Lot: 11710-52  
Prod.: Vet-Con Norm  
ID: 1171052+  
Exp.: 26-Mar-2018

